

We've got you covered!

The widest range of myopia management contact lenses



Daily Disposable

Quarterly Replacement



Night Lenses



www.positiveimpact.co.uk

Telephone: 01424 851375 | Email: info@positiveimpact.co.uk

positive
impact



Daily Disposable

NaturalVue® Enhanced 1Day Multifocal

NaturalVue® is a patented EDOF lens that helps take control of paediatric myopia. When fitted to myopic children it demonstrates an 85% decrease in the rate of myopia progression.¹ CE Marked for presbyopia and myopia control.

Material:	Etafilcon A (58% Water)
Parameters:	Base Curve: 8.3 (nominal) Effective Add up to +3.00D Diameter: 14.50mm
Power Range:	+4.00 to -12.25D (0.25D steps)



Night Lenses

Alexa® AR & Alexa® AR Toric

Alexa® AR is a new concept in OrthoK, which provides faster emmetropisation and optimised memory capacity.

Material:	Paragon HDS 100 (FDA approved for overnight wear)
Diameter:	Alexa AR (all): 9.80 to 12.40mm
Power Range:	Alexa AR: up to -10.00D (0.25 steps) Alexa AR Toric: Cyl upto -3.50D (0.25D steps)
Compression:	up to -5.00D: +1.25D -5.25 to -10.00D: +2.00D



Quarterly Replacement

Esencia® & Esencia® Toric

Esencia® is a clinically proven², customised peripheral addition soft contact lens indicated for the correction and control of myopia progression.

Material:	Acofilcon B (50% Water)
Parameters:	Base Curve: 7.80 to 9.20 (0.10mm steps) Diameter: 14.00 and 14.50mm
Power Range:	Esencia: Plano to -20.00D (0.25D steps) Esencia Toric: Cyl Power: -0.75 to -6.00D (0.25D steps) Cyl Axes: 5 to 180 degrees (5° steps)

www.positiveimpact.co.uk

Telephone: **01424 851375** | Email: info@positiveimpact.co.uk

References:

1. Cooper J, O'Connor B, Aller T, et al. Reduction of myopic progression using a multifocal soft contact lens: A retrospective cohort study. Clin Ophthalmol. 2022Jul;16:2145-2155.
2. García-del Valle, Blázquez, Gros-Otero et al. Efficacy and safety of a soft contact lens to control myopia progression. Clin Exp Optom 2020.

positive
impact