

cleadew  
CareSolution

## Rinsing solution for all rigid gas permeable contact lenses.


Introducing cleadew CareSolution, the perfect rinsing companion for all rigid gas permeable contact lenses.

Designed to seamlessly integrate with other Cleadew lens care products, it ensures a hassle-free and reliable lens care routine. This solution promotes optimal eye health with its gentle yet effective formula. Easy to use and suitable for all lens types, **cleadew CareSolution** ensures a straightforward approach to maintaining clean and hygienic contact lenses.



- 1 cleadew CareSolution:  
360mL x 1  
Hydrogen peroxide as  
preservative (0.004%),  
Boric acid

- 2 cleadew CareSolution  
is used as follows:



Pick up the cleaned / disinfected lenses from the lens case. Hold them with your fingers and rinse both sides with the cleadew CareSolution for at least 10 seconds.