

OphTECS cleadew contact lens care solution



Care system for rigid oxygen-permeable contact lenses (optimal for ortho-K)

- Disinfection, cleaning, rinsing, storage
- Proteolytic enzymes and surfactants remove lipid and protein deposits
- Povidone iodine is highly effective against microorganisms and acathamoeba
- Optimal cleaning without rubbing the contact lenses and easy to use, disinfection completed after 4 hours
- Package contains cleaning solution, neutralizing tablets, rinsing solution* and a contact lens case
- Available in:
 - cleadew GP 40ml & cleadew GP 120ml
 - cleadew GP 40ml + *Care Solution 120ml
 - cleadew GP 120ml + *Care Solution 120ml



All-in-one solution for rigid oxygen-permeable contact lenses

- Disinfection, cleaning, rinsing, storage
- Improved wettability thanks to innovative SMD technology
- Greater wearing comfort thanks to the accumulation of hyaluronic acid derivatives (HAD) on the contact lens surface
- Perfect balance between disinfection and safety for the eye
- Package includes all-in-one solution and a contact lens case
- Available in 60ml & 360ml



Care system for scleral contact lenses

- Disinfection, cleaning, rinsing, storage
- Proteolytic enzymes and surfactants remove lipid and protein deposits
- Povidone iodine is highly effective against microorganisms and acanthamoeba
- Optimal cleaning without having to rub the contact lenses
- Easy to use, disinfection completed after 4 hours
- Package contains cleaning solution, neutralizing tablets, rinsing solution* and a large format contact lens case
- Available in:
 - cleadew SL 100ml & cleadew SL 300ml
 - cleadew SL 100ml + *Care Solution 120ml
 - cleadew SL 300ml + *Care Solution 360ml



Rinsing and insertion solution in ampoules for scleral contact lenses

- Specially designed for rinsing and inserting scleral lenses
- Added sodium hyaluronate increases wearing comfort and enables a longer wearing time
- Thanks to the hypromellose, the contact lens is easier to insert and the risk of air bubbles is reduced
- Preservative-free
- Package contains 30 ampoules of 8ml

OphTECS cleadew contact lens care solution



Care system for soft contact lenses

- Disinfection, cleaning, rinsing, storage
- Proteolytic enzymes and surfactants remove fat and protein deposits
- Povidone iodine is highly effective against microorganisms and acanthamoeba
- Easy to use without rubbing lens, disinfects within 4 hours
- Package contains cleaning and rinsing solution, disinfection and neutralization tablets and a contact lens case
- Available in 128ml & 385ml



All-in-one solution for soft contact lenses

- Disinfection, cleaning, rinsing, storage
- World's first all-in-one contact lens solution with a combination of PHMB and alexidine
- Highly effective against microorganisms and acanthamoeba, reduces the risk of eye infections
- Improved wettability thanks to innovative SMD technology
- Increased wearing comfort thanks to the accumulation of hyaluronic acid derivatives (HAD) on the contact lens surface
- Package includes all-in-one solution and a contact lens case
- Available in 60ml & 360ml



Rinsing solution for all rigid oxygen-permeable contact lenses

- The perfect rinsing solution for hard oxygen permeable contact lenses.
- Designed to work seamlessly with other cleadew lens care products and provides a hassle-free and reliable lens care routine.
- Promotes optimal eye health with its gentle yet effective formula.
- Cleadew CareSolution provides a simple approach to achieve and maintain clean and hygienic contact lenses.
- Available in 360ml